

There are three sizes of racks designed to hold blood bag cassettes. The two vertical racks will hold 11 and 20 cassettes respectively. They are part number 11799031 and 11802121. The third rack will hold 11 bag cassettes horizontally. This part number is 11808961. Call Customer Service for pricing and availability.

You Wanted To Know....

- Q. Is there anyway of changing the 60 minute default time for the LN2 supply without using a computer interface?
- R. No. But by connecting with a PC this parameter can be set anywhere from 1 minute to 360 minutes.
- Q. Is it better to ship units with weight on casters or weight on bottom head?
- R. It is safer to ship them with the weight on the bottom head than on the casters. The impact loads on the casters from riding in a road trailer are probably higher than the load rating on the casters. This is why the heavier models are crated on blocked skids.
- Q. Approximately what year did MVE go to the 1.6 version of the e-prom chip?
- R. The following timeline is close to when different versions of the e-prom chip took effect:
- | | |
|------------------------|--------------------|
| TEC 2000 | |
| Prior and up to 1995- | Lakeshore |
| 1/96 up to 4/1/98- | MDC 1.5 version |
| 4/1/98 up to 1/1/2000- | MDC 1.6 version |
| 1/1/2000 to 1/1/2001- | MDC 1.7 version |
| 1/1/2001 to 4/1/2003- | MDC 1.71 version |
| 4/1/2003 to present- | MDC newest version |

Note to remember

Note to anyone who works with LN2 or LO2. The higher the pressure, the warmer the liquid. The coldest liquid is that which has a low saturation pressure.

Feel free to call Technical service @ 866 819-5897 for any questions concerning Bio-Medical dears or equipment. Contact Customer service for pricing and inventory questions.

Upcoming Show Dates:

ASH (American Society of Hematology)
San Diego CA
Dec 5-9, 2003

Bio-Medical Customer and Technical Service

| | |
|-------------------|---|
| Customer Service | 888 683-2796 toll free / 952 882-5000 Burnsville receptionists 800 232-9683 fax |
| Technical Service | 952 641-6115 direct line 866 819-5897 toll free 612 382-6678 cell 800-232-9683 fax |

For copies of past Tech Tips or for more information on maintaining your nitrogen storage dewars please contact Jim Bachman at (952) 641-6115, Fax (800) 232-9683.

