



# Chart Download Tool Information

	Old Tool	ChartConnect 3000
TEC 3000 Compatibility	Limited	Fully Compatible
TEC 2000 Compatibility	Limited	Fully Compatible
Network Compatible	No	Yes
Automatic Controller Detection	No	Yes
Com Port Status Display	No	Yes
Output File Format	Text File	CSV file (Excel Compatible)
Selectable Download Range	No	Yes
Max Fill Time Status Change	No	Yes
Log Period Status/Change	Yes	Yes
Unit Setting Display	No	Yes
Level Settings Display/Change	No	Yes
Temperature Settings Display/Change	No	Yes
Advanced Time Settings (troubleshooting & Optimization)	No	Yes
Transfer Speed (records/second)	2.0	5.8

Chart Industries, Inc.  
MVE, BioMedical Group  
2200 Industrial Drive, Suite 500  
Ball Ground, GA 30107  
Ph 770-721-7759 • Toll Free 1-800-482-2473  
Fax 770-721-7758



[www.chartbiomed.com](http://www.chartbiomed.com)

Chart BioMedical, Ltd.  
Unit 2, Maxdata Centre  
Downmill Road, Bracknell  
Berks RG12 1QS, United Kingdom  
Ph +44(0) 1344 403100  
Fax +44(0) 1344 429224