



## MVE TEC 2000 vs. TEC 3000

	TEC 2000	TEC 3000
Alarm Outputs	No Discrete Alarm Indicators	Discrete Alarm Indicators Temp A High Low Level High Level Low Battery
Serial Communications	One RS-485 Serial Port	Two RS-485 Serial Ports
User Interface	2 x 16 LCD Display 20 Membrane Switches (Keys)	4 x 20 LCD Display 8 Membrane Switches (Keys)
Data Storage	360 Records	30,000 Records
Security	1 Password	10 Programmable Passwords 4 Possible Security Levels
Timed Filling	N/A	Programmable 1-28 Days
One Fill All Fill	Simultaneous Only	Simultaneous or Sequential
Sequential Fillings	N/A	Capable To Fill Up To 100 Units
Languages	English	English Spanish French German Italian
Bypass Sensor Calibrator	No	Yes
Sensor Calibration Alarm	No	Yes Temp A Temp B Bypass Temp
Usage Warning	No	Yes, When usage doubles from one calculation to the next
Usage Alarm	No	Yes, When usage increases by 5x from one calculation to the next