



MVE Tech Tips



January, 2007


Lab 20 Discharge Device Kit

The lab 20 discharge device has undergone an engineering change. We have an extended stopper available that will replace the existing ones in order to increase efficiency of the device. An upgrade kit for older discharge devices has been created (P/N 13490811). Detailed installation instructions will be included with the upgrade kit.

The new discharge device is part number 13464867, and the old number, 10667802, has been obsolete. Since the Lab 10 discharge device used the same rubber stopper, a new part number was also created for the Lab 10 discharge device. The new number is 13467208; this part number applies to Lab 10 dewars with a 2" neck opening. An updated manual is currently being created. Remember that the warranty on all parts is 90 days.

Clear 2" Boxes

We have begun stocking clear boxes for use with our freezers and some of our aluminum dewars. These boxes are 2" tall square boxes and have a capacity of 100 internally threaded vials. The list price on these is \$10.15. Please use part number 9719649CLR when ordering.



Replacement lids for freezers

As a result of low volume production, manufacturing costs on many of our service lids fluctuate rapidly. Because of this, it is not feasible for us to keep up to date pricing information in our inventory system. To ease the purchasing process we recommend that you call in to receive a quote on price and delivery time before placing an order. This allows customer service to obtain updated cost and deliver times from purchasing to insure an accurate quote. To find applicable part numbers for replacement lids, please see August 2006 Tech Tips.

Customer/Technical Info

For more information please contact Technical Service:

For ordering information contact:

USA & Canada: 800 482 2473
Jennifer Brewer, Wendy Cooke, &
Kari Cockrell

International: 800 482-2473
Nicole Curtis & Heidi Parker

Europe: +44 1344 403100

Asia/Australia: +61 297 494333

Technical Service:

Jim Bachman: 952 758-8411

Mark Baglia: 770 257-1269

w w w . c h a r t b i o m e d . c o m