



MVE Tech Tips



June, 2007

Part number changes

In May of this year, several part numbers across the AI product line were changed. Please be aware of these changes:

Old P/N	Old Description	New P/N	New Description
9723199	canister sq 2" box	11554720	canister x-large box

As a result of the change above, all of the AI units using this part number have been assigned new part numbers.

10811464	Fnl XC47/11-6 w/2"sq canisters	13346449	Fnl XC47/11-6 w/XL sq canisters
13101839	XC34/18 "Duke"	13101839	XC 34/18 Plus w/weld canisters

Fuzzy cork/cover systems

In response to customer demand, "fuzzy" cork and cover assemblies are now available. The cork and cover assemblies below are compatible with cryoshipper units.

11924891	Service C/C with fleece only
13603708	Service C/C with fleece and datalogger attached
13383063	Final CryoShipper NMDP, this is a Cryoshipper with fleece cork and data logger attached

Blood Bag Size and Descriptions

There has been some confusion in the past about blood bag descriptions and correlating sizes. Please see the chart below detailing cassettes and frames provided by Chart.

Blood Bag Sizes and Descriptions		
Bag size	Frame size	Part Number
50 ml Baxter 4R9951		
cassette		10823393
frame	6 place	10823385
	7 place	11831951
	8 place	11831969
250 ml Baxter 4R9953 (also Fenwal 4R5461)		
cassette		10580694
frame	3 place	10580758
	4 place	10580740
	5 place	10968111
	6 place	11550905
500 ml Baxter 4R9955 (also Fenwal 4R5462)		
cassette		10580707
frame	3 place	10580774
	4 place	10580766
	5 place	10968437
750 ml Baxter 4R9957		
cassette		11376281
frame	4 place	11376302
1000 ml Baxter 4R9959		
cassette		11875929
frame	4 place	11875953
200 ml Gambro DF-200		
cassette		10580715
frame	3 place	10580791
	4 place	10580782
	5 place	10969149
700 ml Gambro DF-700		
cassette		10580723
frame	3 place	10580811
	4 place	10580803
	5 place	10823115
1200 ml Gambro DF-1200		
cassette		11041617
frame	4 place	11041625

Bag size	Frame size	Part Number
25 ml Medsep		
cassette		11856322
frame	6 place	11877158
	7 place	11044623
	8 place	11856314
500 ml Charteredmed 2030-2		
cassette		10580707
frame	3 place	10580774
	4 place	10580766
	5 place	10968437
750 ml Charteredmed 3030-2		
cassette		10580731
frame	2 place	10580838
	3 place	10580820
500 ml Charteredmed 7005		
cassette		10969077
frame	3 place	10969069
	4 place	10969114
750 ml Charteredmed 7400		
cassette		10969472
frame	2 place	10969481
	3 place	10969464
Heart valve CBS		
rack only	4 place	11011944

Customer/Technical Info

For ordering information contact

Customer Service:

USA & Canada: 800 482-2473
 Europe: +44 1344-403100
 Asia/Australia: +61 297-494333

For technical information contact

Technical Service:

USA & Canada
 Jim Bachman: 952 758-8411
 jim.bachman@chart-ind.com
 Mark Baglia: 770 257-1269
 mark.baglia@chart-ind.com
 Europe
 Kieran Taheny: +44 1344 403100
 kieran.taheny@chart-ind.com

www.chartbiomed.com