



Tech Tips




September 2010

IATA Cryovials

Many varieties of vials are available to our customers today. Chart offers plastic, externally threaded Cryovials as an IATA approved accessories. Each part number is 100 vials.

Externally Threaded Cryovials

| Description | Part Number |
|------------------------------|-------------|
| CRYOVIAL IATA PACKAGE 1.2 ML | 10778810 |
| CRYOVIAL IATA PACKAGE 2ML | 10778828 |
| CRYOVIAL IATA PACKAGE 3ML | 10778836 |
| CRYOVIAL IATA PACKAGE 4ML | 10777438 |



Manual Discharge Device Transfer Hoses

Several filling accessories are available for the manual discharge device. The stainless steel transfer hoses are two feet long, have threads ¼ inch in diameter, and are equipped with a choice of two phase separators.

Stainless Steel Transfer Hose

| Phase Separator Size | Part Number |
|----------------------|-------------|
| Small | 14044143 |
| Medium | 14044151 |



Tec 2000 to TEC 3000 Upgrade Wiring, PN: 13319512

When wiring your stand alone TEC 3000 upgrade kit for MDC TEC 2000 controllers, please be mindful of the battery backup polarity. Also, be sure you are working from the correct instruction revision. The document should read PN 13319601 Rev A.

Customer & Technical Info

For ordering information contact
Customer Service:

Europe: +44 1344 403100
Australia/New Zealand: +61 2 9749 4333
US Toll Free: 800.482.2473
Asia/Pac Rim: 770.721.7759
The Americas: 770.721.7759

For technical information contact
Technical Service:

Europe:
 Phone: +44 1344 403100
 Email: techservice.europe@chart-ind.com
US Toll Free
 Phone: 800.482.2473
 Email: techservice.usa@chart-ind.com
Asia/Pac Rim
 Phone: 770.721.7759
 Email: techservice.usa@chart-ind.com