

With the PDF Logger®, you require only one single device and no specific software or cables to monitor temperature during the global transfer of cryobiological goods. The incredibly easy-to-use “Transit” and “Arrived” buttons make initiating and receiving monitored shipments fast and simple. The PDF Logger can be plugged into the USB of any PC to automatically generate a PDF evaluation report including statistics, alarm status, and temperature graph which can be read anywhere without difficulty.

With the 1-year lifespan, there are no more hassles of recalibration and recertification. After the datalogger has been in service for 1 year, simply discard the logger and purchase a new one for about the same cost of recalibration and battery replacement. Your PDF Logger is always under warranty!



Specifications

Dimensions L, W, H in (mm)	3.7 x 1.6 x 0.5 (95 x 40 x 12)
Weight oz (g)	1.4 (40)
Temperature Range	-200°C to +200°C
Accuracy	± 0.5°C [-200°C, -10.1°C]
Resolution	0.1°C
Display	Multifunction LCD, 23.5 x 23.5 mm
Case	ABS plastic
Battery	3.6V
Battery Life	400 days
Sensor	Pt100 3-wire class 0.5
Memory	16,000 data points
Interface	USB – PC Universal Serial Bus
Software Compatibility	Windows 7, VISTA, XP, Win2000, Win 98/ME
Configuration	LIBERO Configuration Utility – Free of charge on www.pdf-logger.com
Evaluation Report	Built-in PDF file generator that automatically establishes an evaluation file including graph. Complies to the ISO Standard 19005-1 Document Management for the long-term preservation of electronic documents

Highlighted Features

- Displays temperature and alarm status
- Records and stores up to 16,000 data points
- High and Low Alarm Settings
- No recalibrations
- 1 year lifespan – after 1 year, simply discard and purchase new logger
- Battery life remaining days displayed
- No additional software or cables needed!
- Automatically generates PDF /A evaluation report including graph
- CE marked
- ISO 19005-1
- **Compliant: FDA 21 CFR Part II, GMP, GLP**

1 Year Warranty

Chart BioMedical
 2200 Airport Industrial Dr., Ste. 500
 Ball Ground, GA 30107
 Ph 770-721-7759 • Toll Free 1-800-482-2473
 Fax 770-721-7758
customerservice.usa@chartindustries.com

By ELPRG



www.chartbiomed.com

Chart BioMedical Ltd.
 Ph +44(0) 1344 403100 • Fax +44(0) 1344 429224
customerservice.europe@chartindustries.com
Chart Japan Co., Ltd.
 Ph + (03) 5776 2670 • Fax + (03) 5776 2676
customerservice.japan@chartindustries.com
Chart Australia Pty Ltd.
 Ph + 61 2 974 94 333 • Fax +61 2 974 94 666
customerservice.australia@chartindustries.com